Fitting procedure using Pylorenzmie

1 Fitting procedure using Pylorenzmie

In order to fit an hologram, I used the pylorenzmie model which provides a set of python classes in order to analyse holographic microscopy data.

Pylorenzmie can be download on the David Grier's github repository: https://github.com/davidgrier/pylorenzmie.

What I actually get from the experiments are mp4 movies, in order to analyze them easily, I constructed a wrapper around the pylorenzmie module which can be found on my repository: https://github.com/eXpensia/wraplorenzmie.

This wrapper permits to do the following pipeline:

- Directly load the movies
- Compute the back ground.
- Use the first image in order to get the pre guesses
- Fit the 10 000 first images to determine precisely the radius and index of a particle.
- Use the later information in order to fit the whole movie (and save the data in the same time)

One that done, the trajectory be analyzed separately.

```
import wraplorenzmie.utilities.utilities as utilities
import wraplorenzmie.fits.fit as fit
import imageio
# For Plotting.
import matplotlib.pyplot as plt
import seaborn as sns
import numpy as np
#sns.set(style='white', font_scale=2)
%matplotlib inline
import matplotlib as mpl
```

```
mpl.rcParams["figure.dpi"] = 200
from matplotlib import rc
rc('font', family='serif')
rc('text', usetex=True)
rc('xtick', labelsize='x-small')
rc('ytick', labelsize='x-small')

def cm2inch(value):
    return value/2.54
```


No module named 'pylorenzmie.fitting.cython.cminimizers'

```
[2]: #We load the movie
vid = utilities.

→video_reader("Basler_acA1920-155um__22392621__20200527_162231224.mp4")
```

```
[4]: # We take a look at the first image of the movie
image = vid.get_image(1)
plt.imshow(image,cmap="gray")
```

[4]: <matplotlib.image.AxesImage at 0x1a70b337be0>



```
[5]: # set the background image (it can also be computed using vid.

→get_background method)

vid.number = 125000


vid.background = np.array(imageio.imread("background.tiff"))

#vid.background = vid.get_background(n=50) # n is the number of image to

→use to compute the background

plt.imshow(vid.background,cmap="gray")
```

[5]: <matplotlib.image.AxesImage at 0x1a70bc56490>




- [6]: imageio.imwrite("background.tiff", vid.background) # We save the background → for possible later use.
- [7]: # the normalized image, we can see that their is some movement in the background.

 # This could be avoided by computed the background as a function of the time, if the particle diffuses enough.


 normed_image = utilities.normalize(image, vid.background)

 plt.imshow(normed_image, cmap="gray")

 normed_image = normed_image



[8]: # We found the possition of the particle feature = utilities.center_find(image)[0] utilities.plot_bounding(normed_image,feature)




1.1 Fitting the first image

We fit the first image in order to get the preguess. We first start by croping the hologram.

```
[9]: xc, yc, w, h = feature[0]
x_center = xc
y_center = yc
h=200
im_c = fit.crop(image, int(xc), int(yc), int(h))
bk_c = fit.crop(vid.background, int(xc), int(yc), int(h))
cropped = utilities.normalize(im_c,bk_c, dark_count = np.min(im_c))
cropped = cropped / np.mean(cropped)
plt.imshow(cropped,cmap = "gray")
```

[9]: <matplotlib.image.AxesImage at 0x1a71d7e6d00>




```
[12]: zo = result.result["x"][2]*0.0513
      print(result.result["x"][2]*0.0513)
      print(result.redchi)
      print(result.result["x"])
     11.427616273713154
     7.2459765196825305
     [101.23514587 103.00299474 222.76055114
                                                 1.5310255
                                                               1.58239091
        1.00198476]
     We can plot the result to see if the fit worked properly, and, for a more quantitative com-
```

```
parison we can compute the radial intensity profile of both hologram and compare them.
[13]: center = np.array(np.shape(fitter.image))
[14]: radial_exp, err = radial_profile(fitter.image)
      theo_exp, err = radial_profile(fitter.fitter.model.hologram().
       →reshape(fitter.shape))
      # computing first the holgram using the fit resutlt
[15]: fit_data = {}
      radius_radial = np.arange(len(radial_exp)) * 0.0513
      plt.figure(figsize = (15,15))
      fig = plt.figure(figsize=(cm2inch(8.6),1.65*cm2inch(8.6)))
      fig.subplots_adjust(left=0.14, bottom=.12, right=.99, top=.98)
      plt.subplot(2,2,1)
      plt.imshow(fitter.image, cmap = "gray")
      #plt.title('subplot(2,2,1)')
      fit_data["exp_image"] = fitter.image
      plt.subplot(2,2,2)
      plt.imshow(fitter.fitter.model.hologram().reshape(fitter.shape), cmap = ___

¬"gray")
      frame1 = plt.gca()
      frame1.axes.yaxis.set_ticklabels([])
      fit_data["th_image"] = fitter.fitter.model.hologram().reshape(fitter.
       ⇒shape)
      #plt.title('subplot(2,2,2)')
      plt.subplot(2,2,(3,4))
```

<Figure size 3000x3000 with 0 Axes>



```
[16]: fitter.fit_video(vid = vid, savefile="find_nrfit_result_dur_27052020_n_r_fix_0p0513_wav532.

dat",xc = x ,yc= y, h = 200, n_end=10000,method = "lm")
```

100% 9999/9999 [12:39<00:00, 13.17it/s]

```
[17]: # Since the measurement or not saved into the ram we need to load it
n_r = np.fromfile('find_nr_exame.dat', dtype=np.float64)
n_r = n_r.reshape(len(n_r)//10,10)
r = n_r[:,3]
n = n_r[:,4]
```

1.2 Fitting the n, r distributiton using a KDE estimator


To find the most probable couple of r/n we use a kde estimator using seaborn

```
[18]: import numpy as np
     import scipy.stats as st
     import matplotlib.ticker as ticker
     data = np.random.multivariate_normal((0, 0), [[0.8, 0.05], [0.05, 0.7]],__
      →100)
     x = r[(r>1.5) & (r<1.555)]
     y = n[(r>1.5) & (r<1.555)]
     xmin, xmax = np.min(x), np.max(x)
     ymin, ymax = np.min(y), np.max(y)
     # Peform the kernel density estimate
     xx, yy = np.mgrid[xmin:xmax:100j, ymin:ymax:100j]
     positions = np.vstack([xx.ravel(), yy.ravel()])
     values = np.vstack([x, y])
     kernel = st.gaussian_kde(values)
     f = np.reshape(kernel(positions).T, xx.shape)
     f = f/np.max(f)
```

```
[19]: np.round(np.max(f))
```

[19]: 1.0

```
fig = plt.figure()
fig.subplots_adjust(left=0.16, bottom=.20, right=.99, top=.99)
ax = fig.gca()
#ax.set_xlim(1.505, 1.53)
#ax.set_ylim(1.575, 1.6)
# Contourf plot
cfset = ax.contourf(xx, yy, f, cmap='Blues')
## Or kernel density estimate plot instead of the contourf plot
#ax.imshow(np.rot90(f), cmap='Blues', extent=[xmin, xmax, ymin, ymax])
# Contour plot
cset = ax.contour(xx, yy, f, colors='k', levels=6)
```



```
[21]: print(" n determined with : mu={0}, sigma={1}".format(np.mean(yy[np.

→where(f > 0.1)]), np.std(yy[np.where(f > 0.1)]))

print(" r determined with : mu={0}, sigma={1}".format(np.mean(xx[np.

→where(f > 0.1)]), np.std(xx[np.where(f > 0.1)])))
```

n determined with : mu=1.5851200393768743, sigma=0.003267685282504072 r determined with : mu=1.5181266656310368, sigma=0.00682411690457934

```
[22]: (mu_n, mu_r) = np.mean(yy[np.where(f > 0.1)]), np.mean(xx[np.where(f > 0.1)])
```

1.3 Fitting the whole movie


Now that the measurement of n and r is one we can move on the measurement of the whole trajectory by simply using fitter.fit_video. For demonstration purposes, I only fit here at $\simeq 22$ image per seconds, if can goes up to at least 60 with recent GPU.

100% 9999/9999 [07:24<00:00, 22.47it/s]

1.4 Plot the trajectory

```
[25]: plt.plot(np.arange(len(z))/60, x)
   plt.plot(np.arange(len(z))/60, y)
   plt.ylabel("position $\mathrm{\mu m}$")
   plt.xlabel("time (s)")
```

```
[25]: Text(0.5, 0, 'time (s)')
```



```
[26]: plt.plot(np.arange(len(z))/60, z)
   plt.ylabel("height $\mathrm{\mu m}$")
   plt.xlabel("time (s)")
```

[26]: Text(0.5, 0, 'time (s)')

